



[Ophthalmic Equipment](#) » 
 [Ophthalmic Supply](#) » 
 [Optical Equipment](#) » 
 [Optical Supply](#) »

## Categories

[Ophthalmic Equipment](#)  
[Ophthalmic Supply](#)  
[Optical Equipment](#)  
[Optical Supply](#)

## Search

Go!

Advanced search

Home > Ophthalmic Equipment > Slit Lamps



## Ezer ESL-2600 Slit Lamp



Price: \$3,740.00

Product ID : 100-100-060

Weight: 48.40 lbs

0 review(s)

Like 0

Save

Quantity

1

Add to Cart



Ezer ESL-2600 Slit Lamp



Ezer ESL-2600 Slit Lamp



Ezer ESL-2600 Slit Lamp



Ezer ESL-2600 Slit Lamp

## Description

The Ezer ESL-2600 Slit Lamp is a traditional style slit lamp with a vertical illumination tower and has standard 10X and 16X magnifications. Sleek design provides increased optical efficiency, light transmission, illumination and depth of focus.

The ESL-2600 Slit Lamp presents a stereoscopic optical system by adding fusing images a focal point.

The ESL-2600 brings a two position lever that is located under the optical head that makes it easy to interchange the magnification levels.

A pair of eyepieces are optional in order to change the magnification levels in cases where increased magnification is essential.

The ESL-2600 Slit Lamp is capable of magnifications of 10x and 16x.

### Joystick

Incline it to move the instrument slightly on the horizontal surface and rotate it to adjust the elevation of the microscope.

### Slit width control knob

The slit width is adjustable from 0 to 9 mm.

### Fixation target

An illuminated fixed spot for the patient to look at.

### Base Locking Screw

The base will be locked when this screw is fastened.

### Position Lever

Apertur and slit height control knob  
Rotate this knob to adjust the spot and slit height. Swing the knob horizontally to revolve the slit.

**&#149; Horizontal mark**

When the horizontal center of the patient's eye is aligned with this mark, the elevation of the microscope. Controlled by joystick is also on its center position.

• **There are four folters for selection**

Heat-absorption, grey, red-free and blue.

• **Slit height and aperture display window.**

• **Lamp cap**

• **Chin-rest elevation adjustment knob**

Rotate the knob to adjust the elevation of the chin-rest.

• **Protection cap**

• **Microscope and illumination arm couple bolt**

Fasten this bolt and the illumination arm and the microscope arm can be Rotated together.

The two position lever located under the optical head, makes it easy to interchange the magnification levels.

**Optional accesories**

(Additional purchase)

**Applanation tonometer**

This ESL-2600 slit lamp can be equipped with TN-150, Model R-900 or Model T-900 applanation tonometer to measure the intraocular pressure.

**Product Reviews**

Login or Register to write the first review.

**Products You May Like**



iBEX 5 Step LED Slit Lamp

\$4,995.00



iBEX 2-Step LED Slit Lamp

\$2,995.00



Ezer ESL-5200 Slit Lamp

\$4,240.00



Ezer ESL-700 Hand Held Slit Lamp

\$3,740.00